

Fusion is the merging of product formats into a single innovative portable product. As a stable unit, it can smoothly roll over the text and magnify. Comfortable reading is achieved with the tiltable display. Its unique scroll zoom camera is easily removed so you can view inaccessible text for handwriting. Fusion is your solution for magnification on the move.

- Modes:** Full colour, enhanced positive and negative reading modes and also a range of semi colours. These are controlled by ergonomic buttons on the mouse. Preferred reading modes are stored for ease of use.
- Magnification:** 7" High Brightness display with zoom range from 6x to 16x.
- Overview:** Unique capability to assist with handwriting and distance viewing.
- Scroll zoom:** Unique, easy to use scroll zoom feature that facilitates easy selection of the required magnification.
- Rollers** on the base unit allow easy side to side movement of the unit when reading text.
- The unit folds conveniently into an overall height of just 31 mm and weighs only 1.1 kg.





Compact folding design



Removeable camera



Built-in handwriting feature

## Everything you need is included:

- Fusion device with integrated screen and mouse camera
- Protective carry case
- Internal lithium ion rechargeable batteries
- Mains power cord
- Power supply/charger
- User manual

## Specification

<b>Art. Nr.:</b>	<b>165630</b>
<b>Screen:</b>	7 inches (180 mm) TFT LCD
<b>Magnification:</b>	6 – 16x
<b>Modes:</b>	Full colour. Enhanced black on white. Enhanced white on black. Semi Colours: Broad range of text and background colours.
<b>Weight:</b>	1.1 kg (2 lb 4 oz)
<b>Unit Size:</b>	174 x 182 x 31 mm (7 x 7 x 2 inches) 200 mm (8 inches) high in use
<b>Power:</b>	Worldwide power supply and charger. Typically three hours use from the integrated rechargeable batteries. Charge time with supplied fast charger: 4 hours. Charging battery does not affect usage.

 **ESCHENBACH**

Eschenbach Optik GmbH  
Schopenhauerstraße 10  
90409 Nürnberg  
Telefon +49 (0) 911 3600 -0  
Telefax +49 (0) 911 3600 -358  
Internet [www.eschenbach-optik.com](http://www.eschenbach-optik.com)  
E-Mail [mail@eschenbach-optik.com](mailto:mail@eschenbach-optik.com)