


The new **Prisma AF** introduces a new level of simplicity and ease of use, with the integration of full autofocus technology, to the Prisma product. Now you can enjoy a broad range of magnification levels and a premium image quality with the touch of a button. Simple to connect and intuitive to use, its discrete design and light weight (less than 1.3 Kg), means it can be easily moved or located in a position to best suit your needs. The advanced image processing from Ash enhances the image quality to deliver a well-defined and smoothly magnified image in a variety of modes. Having premium image quality decreases your eye fatigue, so now you can read more quickly and for longer periods of time.



Prisma AF



Reading and writing are so easily achieved with this simple and robust design. Hobbies that had become too difficult to maintain can once more be enjoyed with the **Prisma AF**.



Specifications

Art. No.:	165421
Display:	Uses any standard PAL or NTSC television set with a composite input
Magnification:	2.5 x to 30 x on a 17" TV screen Magnification is dependent on screen size
Modes:	Full-colour Enhanced Black on White, enhanced White on Black
Weight:	1.3 kg
Dimensions:	140 x 82 x 23 mm
Power:	90 - 230Vac with adapter provided
Illumination:	LED illumination The direction of illumination can be varied to suit the user

 **ESCHENBACH**

Eschenbach Optik GmbH
Schopenhauerstraße 10
90409 Nürnberg
Telefon +49 (0) 911 3600 -0
Telefax +49 (0) 911 3600 -358
Internet www.eschenbach-optik.de
E-Mail mail@eschenbach-optik.de