

Video magnifier

Prisma is a stylish electronic reading device which can be connected directly to a standard television and allows you to view text in enhanced reading modes* or photos in full colour on the television screen*. **Prisma** is a versatile, portable product which folds into a compact, transportable configuration weighing just over 1 kg. The **Prisma** can be used to carry out a wide variety of tasks such as needlework and drawing.

- *Modes:** 4-colour photo mode and enhanced positive and negative enhanced reading modes.
- Magnification:** Controlled by adjusting the height of the camera arm as needed. Capable of up to 35 x magnification on a 14" TV screen.
- Focus control:** An easy to adjust focus knob refocuses the image.
- Brightness control:** This varies the brightness in negative mode.
- Illumination:** LED illumination can be switched on in low ambient light conditions.
- Writing mat:** Makes writing with Prisma even easier.





Compact folding design for transport



Variable magnification



Everything you need is included:

- Prisma device and writing pad
- Power supply with connection cable (230 V, 50 Hz)
- User manual
- Protective carrying case

Technical Information

Item No.:	165420
Magnification:	14" screen: 3.8x - 35x 21" screen: 5x - 52x
Modes:	4-colour in photo mode Enhanced black/white Enhanced white/black
Power requirement:	230 V, 50 Hz
Weight:	1.2 kg
Dimensions:	240 x 410 mm
Display required:	Standard PAL television with scart or composite video input (FBAS signal)

ESCHENBACH

Eschenbach Optik GmbH
Schopenhauerstraße 10
90409 Nürnberg
Telefon +49 (0) 911 3600 -0
Telefax +49 (0) 911 3600 -358
Internet www.eschenbach-optik.de
E-Mail mail@eschenbach-optik.de