

OPTi Lite combines the freedom of your laptop with our popular Prisma magnifier to create a versatile mini desktop-style solution. Ideal for both reading and writing, **OPTi Lite** has a wide range of magnification, fine focus control and a 10 hour battery life that won't deplete your laptop. With its stylish carrycase and super light weight, **OPTi Lite** allows you to work wherever and whenever you choose.

- Fully compatible with AI Squared and Dolphin screen magnification software*.
- Light weight, collapsable design that stows in its carrycase.
- Front buttons for fast access to display modes.
- Adjustable magnification and focus.
- Hi-Speed USB 2.0 interface with rock solid, high-resolution, progressive scan image



Supplied with our third generation of OPTi software:

- Unique toolbar interface
- Fully adjustable image window size
- Full colour and enhanced display modes
- Image capture and save facility
- Full control of video display attributes
- Aspect ratio control

Everything you need is included:

- Camera unit with integrated table base
- USB lead for connection to the computer
- OPTi software – powerful and customizable, to display the image on the computer
- Mains power cord
- Power supply/charger
- Carrycase

* Minimum version requirements for compatibility:

- OPTi version 1.12
- ZoomText version 9.0.4
- or
- Lunar/Supernova version 7.0.2

Specification

Art. No.:	165423
Description:	USB 2.0 camera with base suitable for reading and writing
Magnification:	Fully adjustable: 4 x - 40 x on a typical 15 inch screen
Modes:	Full colour Complete range of enhanced colour options Contrast enhancement black/white and inverse
Power:	Worldwide power supply via universal power-supply unit Built-in rechargeable Li-ion battery 10 hours of run time from a charge Charge time with supplied fast charger: 3-4 hours
Weight:	1200 grams (2 lbs 10 oz)
Unit size:	410 x 240 mm base. 280 mm high, collapsing to 70 mm for storage.
System requirements:	Operating System: Windows XP Home or Pro with Service Pack 2/Windows Vista Hardware: 480 MB of RAM (minimum) 640 MB of video RAM minimum; 128 MB or greater recommended

 **ESCHENBACH**

Eschenbach Optik GmbH
Schopenhauerstraße 10
90409 Nürnberg
Telefon +49 (0) 911 3600 -0
Telefax +49 (0) 911 3600 -358
Internet www.eschenbach-optik.de
E-Mail mail@eschenbach-optik.de